

# STAT6 Antibody



PACO01552

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IP

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug

**Protein Background:****Gene ID:**

STAT6

**Uniprot**

P42226

**Synonyms:**

STAT6; Signal transducer and activator of transcription 6; IL-4 Stat

**Immunogen:**

Synthesized peptide derived from human Stat6 around the non-phosphorylation site of T645.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of HeLa cells using Stat6 Polyclonal Antibody.

