

STMN3 Antibody



PACO01553

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

STMN3

Uniprot

Q9NZ72

Synonyms:

STMN3; SCLIP; Stathmin-3; SCG10-like protein

Immunogen:

Synthesized peptide derived from the Internal region of human Stathmin-3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 293 cells using Stathmin-3 Polyclonal Antibody.

