

SYT13 Antibody



PACO01565

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

SYT13

Uniprot

Q7L8C5

Synonyms:

SYT13; KIAA1427; Synaptotagmin-13; Synaptotagmin XIII; SytXIII

Immunogen:

Synthesized peptide derived from the Internal region of human Synaptotagmin XIII.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HuvEc cells using Synaptotagmin XIII Polyclonal Antibody.

