

TNKS1BP1 Antibody



PACO01578

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

TNKS1BP1

Uniprot

Q9C0C2

Synonyms:

TNKS1BP1; KIAA1741; TAB182; 182 kDa tankyrase-1-binding protein

Immunogen:

Synthesized peptide derived from the C-terminal region of human TAB182.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 293 cells using TAB182 Polyclonal Antibody.

