

TSPAN8 Antibody



PACO01607

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

TSPAN8

Uniprot

P19075

Synonyms:

TSPAN8; TM4SF3; Tetraspanin-8; Tspan-8; Transmembrane 4 superfamily member 3;
Tumor-associated antigen CO-029


Immunogen:

Synthesized peptide derived from the Internal region of human Tetraspanin-8.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using Tetraspanin-8 Polyclonal Antibody.