

# GTF2E2 Antibody



PACO01611

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

GTF2E2

**Uniprot**

P29084

**Synonyms:**

GTF2E2; TF2E2; Transcription initiation factor IIE subunit beta; TFIIE-beta; General transcription factor IIE subunit 2

**Immunogen:**

Synthesized peptide derived from the Internal region of human TFIIE-Beta.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of Jurkat cells using TFIIIE-Beta Polyclonal Antibody.

