

TIGD3 Antibody



PACO01626

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

TIGD3

Uniprot

Q6B0B8

Synonyms:

TIGD3; Tigger transposable element-derived protein 3

Immunogen:

Synthesized peptide derived from the C-terminal region of human TIGD3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using TIGD3 Polyclonal Antibody.