

TIMP4 Antibody



PACO01630

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

TIMP4

Uniprot

Q99727

Synonyms:

TIMP4; Metalloproteinase inhibitor 4; Tissue inhibitor of metalloproteinases 4; TIMP-4


Immunogen:

Synthesized peptide derived from the C-terminal region of human TIMP-4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAT-KIDNEY RAT-MUSCLE cells using TIMP-4 Polyclonal Antibody.