

# TLE4 Antibody



PACO01633

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

TLE4

**Uniprot**

Q04727

**Synonyms:**

TLE4; KIAA1261; Transducin-like enhancer protein 4

**Immunogen:**


Synthesized peptide derived from the Internal region of human TLE4.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of RAT-MUSCLE cells using TLE4 Polyclonal Antibody.