

TMEM145 Antibody



PACO01636

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

TMEM145

Uniprot

Q8NBT3

Synonyms:

TMEM145; Transmembrane protein 145

Immunogen:

Synthesized peptide derived from the Internal region of human TMEM145.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HepG2 cells using TMEM145 Polyclonal Antibody.

