

PACO01651

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

TSN

Uniprot

Q15631

Synonyms:

TSN; Translin; Component 3 of promoter of RISC; C3PO

Immunogen:

Synthesized peptide derived from the Internal region of human Translin.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of Jurkat cells using Translin Polyclonal Antibody.

