

TUBA3C/TUBA3E Antibody



PACO01664

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

TUBA3C/TUBA3E

Uniprot

Q13748/Q6PEY2

Synonyms:

TUBA3C; TUBA2; TUBA3D; Tubulin alpha-3C/D chain; Alpha-tubulin 2; Alpha-tubulin 3C/D; Tubulin alpha-2 chain; TUBA3E; Tubulin alpha-3E chain; Alpha-tubulin 3E


Immunogen:

Synthesized peptide derived from the C-terminal region of human Tubulin alpha-3C/D/E.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 cells using Tubulin alpha-3C/D/E Polyclonal Antibody.