

TYK2 Antibody



PACO01666

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

TYK2

Uniprot

P29597

Synonyms:

TYK2; Non-receptor tyrosine-protein kinase TYK2


Immunogen:

Synthesized peptide derived from human Tyk 2 around the non-phosphorylation site of Y1054.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using Tyk 2 Polyclonal Antibody.