

# VAV1 Antibody



PACO01691

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

VAV1

**Uniprot**

P15498

**Synonyms:**

VAV1; VAV; Proto-oncogene vav

**Immunogen:**


Synthesized peptide derived from human Vav1 around the non-phosphorylation site of Y174.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of Jurkat cells using Vav1 Polyclonal Antibody.