

TNFSF15 Antibody



PACO01695

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

TNFSF15

Uniprot

O95150

Synonyms:

TNFSF15; TL1; VEG1; Tumor necrosis factor ligand superfamily member 15; TNF ligand-related molecule 1; Vascular endothelial cell growth inhibitor


Immunogen:

Synthesized peptide derived from the C-terminal region of human VEG1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using VEGI Polyclonal Antibody.