

YEATS2 Antibody



PACO01712

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

YEATS2

Uniprot

Q9ULM3

Synonyms:

YEATS2; KIAA1197; YEATS domain-containing protein 2

Immunogen:

Synthesized peptide derived from the Internal region of human YEATS2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A