

# ZADH2 Antibody



PACO01714

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

ZADH2

**Uniprot**

Q8N4Q0

**Synonyms:**

ZADH2; Zinc-binding alcohol dehydrogenase domain-containing protein 2

**Immunogen:**


Synthesized peptide derived from the Internal region of human ZADH2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of COLO205 cells using ZADH2 Polyclonal Antibody.