

ZNF771 Antibody



PACO01734

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

ZNF771

Uniprot

Q7L3S4

Synonyms:

ZNF771; Zinc finger protein 771; Mesenchymal stem cell protein DSC43


Immunogen:

Synthesized peptide derived from the Internal region of human ZNF771.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using ZNF771 Polyclonal Antibody.