

B3GALT1 Antibody



PACO01738

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

B3GALT1

Uniprot

Q9Y5Z6

Synonyms:

B3GALT1; Beta-1; 3-galactosyltransferase 1; Beta-1,3-GalTase 1; Beta3Gal-T1; Beta3GalT1; UDP-galactose: beta-N-acetyl-glucosamine-beta-1,3-galactosyltransferase 1


Immunogen:

Synthesized peptide derived from the Internal region of human Beta-1,3-Gal-T1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using Beta-1,3-Gal-T1 Polyclonal Antibody.