

B3GALT4 Antibody



PACO01739

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

B3GALT4

Uniprot

O96024

Synonyms:

B3GALT4; GALT4; Beta-1; 3-galactosyltransferase 4; Beta-1,3-GalTase 4; Beta3Gal-T4; Beta3GalT4; GalT4; b3Gal-T4; Gal-T2; Ganglioside galactosyltransferase; UDP-galactose: beta-N-acetyl-galactosamine-beta-1,3-galactosyltransferase


Immunogen:

Synthesized peptide derived from the Internal region of human Beta-1,3-Gal-T4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 cells using Beta-1,3-Gal-T4 Polyclonal Antibody.