

# BAG5 Antibody



PACO01780

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

BAG5

**Uniprot**

Q9UL15

**Synonyms:**

BAG5; KIAA0873; BAG family molecular chaperone regulator 5; BAG-5; Bcl-2-associated athanogene 5

**Immunogen:**


Synthesized peptide derived from the N-terminal region of human BAGE5.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of Jurkat K562 cells using BAGE5 Polyclonal Antibody.