

# RAB26 Antibody



PACO01793

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

RAB26

**Uniprot**

Q9ULW5

**Synonyms:**

RAB26; Ras-related protein Rab-26

**Immunogen:**


Synthesized peptide derived from the C-terminal region of human Rab 26.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of K562 cells using Rab 26 Polyclonal Antibody.