

# STRN Antibody



PACO01810

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

STRN

**Uniprot**

O43815

**Synonyms:**

STRN; Striatin

**Immunogen:**


Synthesized peptide derived from the C-terminal region of human Striatin.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of K562 Jurkat cells using Striatin Polyclonal Antibody.