

CYP2A6/CYP2A7/CYP2A13 Antibody



PACO01811

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

CYP2A6/CYP2A7/CYP2A13

Uniprot

P11509/P20853/Q16696

Synonyms:

CYP2A6; CYP2A3; Cytochrome P450 2A6; 1,4-cineole 2-exo-monooxygenase; CYP2A6; Coumarin 7-hydroxylase; Cytochrome P450 IIA3; Cytochrome P450(I); CYP2A7; Cytochrome P450 2A7; CYP2A7; Cytochrome P450 IIA4; CYP2A13; Cytochrome P450 2A13; CYP2A13


Immunogen:

Synthesized peptide derived from the Internal region of human CYP2A6/7/13.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using CYP2A6/7/13 Polyclonal Antibody.