

CA4 Antibody



PACO01842

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

CA4

Uniprot

P22748

Synonyms:

CA4; Carbonic anhydrase 4; Carbonate dehydratase IV; Carbonic anhydrase IV; CA-IV


Immunogen:

Synthesized peptide derived from the Internal region of human CA IV.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of rat muscle cells using CA IV Polyclonal Antibody.