

BNIP3 Antibody



PACO01868

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

BNIP3

Uniprot

Q12983

Synonyms:

BNIP3; NIP3; BCL2/adenovirus E1B 19 kDa protein-interacting protein 3


Immunogen:

Synthesized peptide derived from the Internal region of human BNIP-3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of mouse eye cells using BNIP-3 Polyclonal Antibody.