

FGF12 Antibody



PACO01875

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

FGF12

Uniprot

P61328

Synonyms:

FGF12; FGF12B; FHF1; Fibroblast growth factor 12; FGF-12; Fibroblast growth factor homologous factor 1; FHF-1; Myocyte-activating factor


Immunogen:

Synthesized peptide derived from the Internal region of human FGF-12.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 22RV-1 BT474 SHSY5Y cells using FGF-12 Polyclonal Antibody.