

FGF18 Antibody



PACO01895

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

FGF18

Uniprot

O76093

Synonyms:

FGF18; Fibroblast growth factor 18; FGF-18; zFGF5


Immunogen:

Synthesized peptide derived from the C-terminal region of human FGF-18.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of SW480 SKOV3 BT474 mouse heart cells using FGF-18 Polyclonal Antibody.