

PACO02062

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

TFRC

Uniprot

P02786

Synonyms:

TFRC; Transferrin receptor protein 1; TR; TfR; TfR1; Trfr; T9; p90; CD71


Immunogen:

Synthesized peptide derived from the N-terminal region of human CD71.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using CD71 Polyclonal Antibody.