

KLF1 Antibody



PACO02106

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

KLF1

Uniprot

Q13351

Synonyms:

KLF1; EKLF; Krueppel-like factor 1; Erythroid krueppel-like transcription factor; EKLF

Immunogen:

Synthesized peptide derived from human EKLF around the non-acetylation site of K274.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A