

# WDR83 Antibody



PACO02121

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

WDR83

**Uniprot**

Q9BRX9

**Synonyms:**

WDR83; MORG1; WD repeat domain-containing protein 83; Mitogen-activated protein kinase organizer 1; MAPK organizer 1

**Immunogen:**


Synthesized peptide derived from the Internal region of human MAPK Organizer 1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of MCF7 cells using MAPK Organizer 1 Polyclonal Antibody.