

RAB5A Antibody



PACO02132

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

RAB5A

Uniprot

P20339

Synonyms:

RAB5A; RAB5; Ras-related protein Rab-5A


Immunogen:

Synthesized peptide derived from the Internal region of human Rab 5A.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Rab 5A Polyclonal Antibody.