

# APOC3 Antibody



PACO02216

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

APOC3

**Uniprot**

P02656

**Synonyms:**

APOC3; Apolipoprotein C-III; Apo-CIII; ApoC-III; Apolipoprotein C3

**Immunogen:**


Synthesized peptide derived from the C-terminal region of human ApoC-III.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of 293T-UV cells using ApoC-III Polyclonal Antibody.