

# DDIAS Antibody



PACO02229

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

DDIAS

**Uniprot**

Q8IXT1

**Synonyms:**

NOXIN; C11orf82; Nitric oxide-inducible gene protein

**Immunogen:**


Synthesized peptide derived from the Internal region of human NOXIN.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of 823 293T AD293 HeLa HepG2 cells using NOXIN Polyclonal Antibody.