

IGFBP1 Antibody



PACO02232

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

IGFBP1

Uniprot

P08833

Synonyms:

IGFBP1; IBP1; Insulin-like growth factor-binding protein 1; IBP-1; IGF-binding protein 1; IGFBP-1; Placental protein 12; PP12

Immunogen:

Synthesized peptide derived from the Internal region of human IGFBP1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of MCF7 cells using IGFBP1 Polyclonal Antibody.

