

TLN1 Antibody



PACO02241

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

TLN1

Uniprot

Q9Y490

Synonyms:

TLN1; KIAA1027; TLN; Talin-1


Immunogen:

Synthesized peptide derived from the C-terminal region of human Talin-1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using Talin-1 Polyclonal Antibody.