

ICAM2 Antibody



PACO02254

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ICAM2

Uniprot

P13598

Synonyms:

ICAM2; Intercellular adhesion molecule 2; ICAM-2; CD102


Immunogen:

Synthesized peptide derived from the Internal region of human ICAM-2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of MCF-7 cells using ICAM-2 Polyclonal Antibody.