

# HTR5B Antibody



PACO02280

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

Htr5b

**Uniprot**

P31387

**Synonyms:**

5-hydroxytryptamine receptor 5B

**Immunogen:**


Synthesized peptide derived from the C-terminal region of human SR-5B.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of COLO cells using SR-5B Polyclonal Antibody.