

CDH3 Antibody



PACO02283

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

CDH3

Uniprot

P22223

Synonyms:

CDH3; CDHP; Cadherin-3; Placental cadherin; P-cadherin


Immunogen:

Synthesized peptide derived from the N-terminal region of human P-cadherin.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using P-cadherin Polyclonal Antibody.