

# CDH11 Antibody



PACO02294

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

CDH11

**Uniprot**

P55287

**Synonyms:**

CDH11; Cadherin-11; OSF-4; Osteoblast cadherin; OB-cadherin

**Immunogen:**


Synthesized peptide derived from the Internal region of human OB-cadherin.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of Jurkat cells using OB-cadherin Polyclonal Antibody.