

VTI1B Antibody



PACO02326

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

VTI1B

Uniprot

Q9UEU0

Synonyms:

VTI1B; VTI1; VTI1L; VTI1L1; VTI2; Vesicle transport through interaction with t-SNAREs homolog 1B; Vesicle transport v-SNARE protein Vti1-like 1; Vti1-rp1


Immunogen:

Synthesized peptide derived from the Internal region of human v-SNARE Vti1p.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAT-HEART cells using v-SNARE Vti1p Polyclonal Antibody.