

VPS26B Antibody



PACO02419

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

VPS26B

Uniprot

Q4G0F5

Synonyms:

VPS26B; Vacuolar protein sorting-associated protein 26B; Vesicle protein sorting 26B


Immunogen:

Synthesized peptide derived from the Internal region of human VPS26B.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of LOVO cells using VPS26B Polyclonal Antibody.