

OR5P2 Antibody



PACO02442

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR5P2

Uniprot

Q8WZ92/Q3MIS8

Synonyms:

Olfactory receptor 5P2


Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 5P2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HUVEC cells using Olfactory receptor 5P2 Polyclonal Antibody.