

OR8I2 Antibody



PACO02445

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR8I2

Uniprot

Q8N0Y5

Synonyms:

OR8I2; Olfactory receptor 8I2; Olfactory receptor OR11-170


Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 8I2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Olfactory receptor 812 Polyclonal Antibody.