

OR2AP1 Antibody



PACO02448

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR2AP1

Uniprot

Q8NGE2

Synonyms:

OR2AP1; OR2AP1P; Olfactory receptor 2AP1; Olfactory receptor OR12-9


Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2AP1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using Olfactory receptor 2AP1 Polyclonal Antibody.