

OR4L1 Antibody



PACO02456

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR4L1

Uniprot

Q8NH43

Synonyms:

OR4L1; OR4L2P; Olfactory receptor 4L1; Olfactory receptor 4L2; Olfactory receptor OR14-28


Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 4L1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using Olfactory receptor 4L1 Polyclonal Antibody.