

OR10H4 Antibody



PACO02466

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR10H4

Uniprot

Q8NGA5

Synonyms:

OR10H4; Olfactory receptor 10H4; Olfactory receptor OR19-28


Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 10H4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAW264.7 cells using Olfactory receptor 10H4 Polyclonal Antibody.