

CNR1 Antibody



PACO02468

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

CNR1

Uniprot

P21554

Synonyms:

CNR1; CNR; Cannabinoid receptor 1; CB-R; CB1; CANN6

Immunogen:

Synthesized peptide derived from the Internal region of human CB1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HT-29 cells using CB1 Polyclonal Antibody.

