

# OR2T11 Antibody



PACO02469

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR2T11

**Uniprot**

Q8NH01

**Synonyms:**

OR2T11; Olfactory receptor 2T11; Olfactory receptor OR1-65

**Immunogen:**

Synthesized peptide derived from the Internal region of human Olfactory receptor 2T11.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of MCF-7 cells using Olfactory receptor 2T11 Polyclonal Antibody.

