

# OR10T2 Antibody



PACO02471

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

OR10T2

**Uniprot**

Q8NGX3

**Synonyms:**

OR10T2; Olfactory receptor 10T2; Olfactory receptor OR1-3

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 10T2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of 293T-UV cells using Olfactory receptor 10T2 Polyclonal Antibody.

